



# ResMed

Actual Size: 5 1/2" (Width)  
x 6 7/8" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

## AirFit™ N20

Nasal mask  
Nasal mask **for Her**



Medical products at home

We provide the official  
document of the manufacturer

[www.multidoctorshop.com](http://www.multidoctorshop.com)

© 2016 ResMed Ltd. All Rights Reserved.  
Specifications may change without notice.  
Product availability may vary across regions.  
1019537/1 2016-07

Rx Only  0123

