

# CAIRE Oxygen Therapy Devices— Treat **all** patients at **any** activity level

	Available Concentrator	Weight	Flow Settings	Oxygen Output	Features
Early Mild	AirSep Focus	1.75 lb	Pulse: 2	332.5 mL/min	<ul style="list-style-type: none"> <li>• Wearable</li> <li>• UltraSense technology*</li> <li>• Simple 1-Touch Operation</li> </ul>
	AirSep FreeStyle	4.4 lbs	Pulse: 1 - 3	500 mL/min	
	AirSep FreeStyle 5	6.2 lbs	Pulse: 1 - 5	1000 mL/min	
Moderate Severe	SeQual eQuinox	14.0 lbs	Pulse: 1 – 9 Cont.: 0.5 – 3.0 LPM	Bolus Size 16 mL – 192 mL	<ul style="list-style-type: none"> <li>• Transportable</li> <li>• 24/7 Capability</li> <li>• autoSAT Technology**</li> <li>• Adjustable Rise Time</li> <li>• Adjustable Sensitivity</li> <li>• Voice Module available on eQuinox</li> </ul>
	SeQual Eclipse	18.4 lbs			

### \*UltraSense® Technology

- -0.054 to -0.1 cm H<sub>2</sub>O sensitivity is the most sensitive trigger of any concentrator
- Prevents skipped breaths and provides maximum patient comfort

### \*\*autoSAT® Technology

- Servo-controlled compressor adjusts speed based on patient's breathing rate
- Keeps saturation levels up as activity increases and ensures fixed bolus size is always delivered with every breath

For more information, please contact our CAIRE Customer Support Team:  
**800.482.2473** [www.chartindustries.com/Respiratory-Healthcare](http://www.chartindustries.com/Respiratory-Healthcare)