

Ambulation & Mobility

fingertip pulse oximeter 7568-SEQ

- Dual-color OLED display
- Indicated for pediatric and adult patients
- Neck strap and batteries included
- Prescription required
- FDA approved
- 1-year warranty



eclipse 3 AC adapter 5941-SEQ



- Includes NEMA power cord
- Intended for use at home, or wherever standard AC power is available
- SeQual branded to help distinguish from E1 and E2 power supplies
- Green LED illuminates when AC Power Supply is supplying 28VDC power

eclipse universal cart 5991-SEQ

- Six-inch wheels for higher ground clearance
- Ideal for all seasons, all terrain
- Easy access to Battery
- Handle slides up and down with different settings for personalized adjustment
- **File HCPCS code E1354**



eclipse wheelchair pack 5220-SEQ



- Attaches to any push-handle wheelchair using adjustable, buckle-style clips
- Eclipse fits easily into the main compartment
- Velcro flap provides easy access to the Eclipse Battery (#7082)



Visit www.sequal.com to view our entire accessory line.
Call 800.482.2473 to order today!



CAIRE
A Chart Industries Company
AirSep® / SeQual®
www.CAIREmedical.com

CAIRE, Inc.
2200 Airport Industrial Way,
Suite 500 Ball Ground, GA 30107
800.482.2473 Telephone
888-932-2473 Fax

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham,
Berkshire, RG41 2PL
Ph +44(0) 1189 367060
Fax +44(0) 0118 9799245

Copyright © 2013 Chart Industries
CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

CAIRE
A Chart Industries Company
AirSep® / SeQual®
www.CAIREmedical.com

MultiDoctorShop
Medical products at home
We provide the official document of the manufacturer
www.multidoctorshop.com

Eclipse Accessories

to support an active lifestyle

Ambulation & Mobility • Travel • Auxiliary & More



Ambulation & Mobility

power cartridge (battery) 7082-SEQ

- Provides up to 5.0 hours of power
- Charges using AC/DC power, or the Eclipse Desktop Charger (#7112)
- Comes standard with each Eclipse
- **File HCPCS code E1356**



desktop charger 7112-SEQ



- Charges the Eclipse Battery (#7082) in less than 3 hours
- Use anywhere in the world, by combining with the Universal Power Plug Set (#5591 sold separately)
- **File HCPCS code E1357**

THE ART OF
OXYGEN

Travel Accessory Kit 5093-SEQ



- Our featured bundle accessory
- Contents include Travel Case (#5010), Battery (#7082), Desktop Charger (#7112) and North American Power Cord
- **File HCPCS codes E1356 (battery) and E1357 (charger)**

Eclipse Protective Cover Item #5052-SEQ

- Reduces external wear and tear
- Protects against spillage and inclement weather
- Secures with convenient zipper closure
- Clear plastic window allows access to control panel



Universal Power Plug Set 5591-SEQ



- Ideal for international travel with the Eclipse
- Europe, UK, New Zealand, North/South America, China, Japan and Taiwan
- Rated at 10A to support the current demand of the Eclipse



Travel Case 5010-SEQ



- Features adjustable Velcro® inner partitions
- Neatly stores Eclipse accessories
- Rubber handle offers a sturdy grip
- Easily slips over the handle of all Eclipse Carts

Eclipse 3 DC Adapter 5942-SEQ

- Intended for use with DC accessory outlets, such as those found in motor vehicles
- Green LED illuminates when DC Power Supply is supplying 26VDC power
- **File HCPCS code E1358**



Eclipse Accessory Bag 7104-SEQ



- Neatly stores Eclipse accessories
- Rubber handle offers a sturdy grip
- Easily slips over the handle of all Eclipse Carts



Eclipse Data Acquisition Tool (EDAT) 5535-SEQ

- Enables firmware upgrades on-site
- Reduces downtime and maintenance costs
- Triages the device for equipment fault vs. user error
- Print and retain historical event logs



Eclipse 2 Plus Firmware Upgrade 5978-SEQ

eclipse² Plus
with *autoSAT*

- Upgrade your Eclipse 2 inventory
- Latest features and technology—same as E3
- Low cost way to stock the latest equipment
- EDAT is required for all upgrades

Humidifier Adapter Kit 7116-SEQ



- Attaches humidifier bottle to Eclipse unit (humidifier bottle not included)