

Effortless
and **effective**



One fit. One solution. One for all.

The perfect balance of sophisticated design and easy functionality, Mirage™ FX for Her is the solution you've been waiting for: performance, ease-of-use and quality all rolled into one.

Benefits for you ...

Sizing is a one-step process, with Mirage FX for Her covering a wide range of facial characteristics in one mask ... so whoever your patients are, Mirage FX is the one for all.

Fitting is just as effortless. Mirage FX is intuitive to fit and adjust, so fitting takes virtually no time at all. With only four parts, including headgear, you'll also save time on managing inventory.

Mirage FX for Her—puts your business first.

... and your patients

The Spring Air™ cushion provides even pressure distribution and comfort for your female patients, incorporating extra softness at the nasal bridge to reduce discomfort commonly felt in this area.

Your patients will also find Mirage FX easy to live with: it's easy to clean, wear and adjust ... and being one of the lightest nasal masks on the market, it won't weigh her down.

Mirage FX for Her—make it your number one mask.



Mirage™ FX for Her
NASAL MASK

www.resmedmiragefx.com

Mirage FX for Her is everything your female patients have been waiting for.

effortless fit


Intuitive adjustments save time and mean just about anyone can achieve a perfect, first-time fit.

simple solution

Mirage FX for Her provides your female patients with performance, ease-of-use and quality all rolled into one.

mask for all

With Mirage FX for Her, you'll fit more people the first time round, saving time, money, and additional patient visits.



New, form-fitting headgear springs to life, ready to fit and wear for easy hair maintenance

SoftEdge™ headgear design is kind to the face, reducing marks and discomfort

Squeeze-tab elbow makes it easy to attach and detach the elbow from the frame

Flexible, winged forehead support allows a dynamic range of fit and adjustment

Lightweight, durable frame is strong enough to stabilise the mask on the face and soft enough to accommodate user movement during sleep

Dual-wall, Spring Air cushion provides even pressure distribution, with extra softness at the nasal bridge, reducing common discomfort in this sensitive area

Diffuse vent design disperses air gently and quietly, away from the bed partner

Headgear loops simplify attachment and detachment of headgear; retain optimal headgear settings

Complete System

	Europe 1	Europe 2	Europe 3
Mirage FX for Her	62132	62133	62139

Europe 1: English (UK), French, German and Italian (Switzerland)

Europe 2: English, Danish, Finnish, Norwegian and Swedish

Europe 3: English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish

Mirage™ FX for Her
NASAL MASK

www.resmedmiragefx.com