



Nebulization simplified

Philips Respironics Family Silver nebulizer compressor system – affordable quality nebulizer therapy

The Philips Respironics Family Silver nebulizer compressor system is a mains powered device ideal for home use

This compact unit is suitable for nebulizing the majority of commonly prescribed aerosol medications and is equipped with a JetBasic nebulizer delivering a quality aerosol.

The Family Silver has a sleek, stylish finish and is supplied with a complete range of accessories including a mouthpiece, nosepiece, adult mask and tubing. In addition a carrying case is provided for convenient storage for patients at home.

Key features

- High flow rate (5 lpm) provides efficient nebulization
- Matched with JetBasic nebulizer for optimum system performance
- 3 year warranty

PHILIPS
RESPIRONICS
sense and simplicity



Simplifying healthcare in the home

Offering high quality performance to help deliver optimum care to patients at home.



Product specification

Weight	1.55 kg
Dimensions (cm)	16.7 x 14.5 x 9.7
Electrical requirements	230V/50Hz/1A (part no.90543) 220V/60 Hz (part no. RDD91099) 220V/60 Hz (part no. 91730) 110V/60Hz (part no. RDD91093)
Maximum nebulizer capacity	6 ml
Noise level	58 dBA
Treatment time (2.5 ml charge)	7 mins
Dynamic flow	5 lpm
Max pressure	2.00 bar
Operating pressure	0.80 bar
Neb rate	0.26 ml/min
MMAD	4.54 µm
Warranty	3 years

The Family Silver has instructions for use in the following languages: English, Français, Italiano, Deutsch, Português, Español, Norsk, Nederlands, Polski, Ελληνικά, Türkçe, العربية, 中文, 日本語, 한국어, Русский, فارسی, Magyar, Čeština, Български, Română, Slovensky and Slovenščina.

Family Silver part numbers

Product description	Part number
Family Silver nebulizer compressor (230V/50 Hz Euro power cord)	90543
Family Silver nebulizer compressor (220V/60 Hz Euro power cord)	RDD91099
Family Silver nebulizer compressor (220V/60Hz Korean power cord)	91730
Family Silver nebulizer compressor (110V/60 Hz Euro power cord)	RDD91093

Family Silver spares and accessories

Product description	Part number
JetBasic nebulizer handset kit containing JetBasic nebulizer, adult mask, pediatric mask, mouthpiece and supply tubing	90326
JetBasic nebulizer handset kit containing JetBasic nebulizer, adult mask and supply tubing	90327
JetBasic nebulizer handset kit containing JetBasic nebulizer, pediatric mask and supply tubing	90328
JetBasic nebulizer handset kit containing JetBasic nebulizer, mouthpiece and supply tubing	90329
JetBasic nebulizer handset only	RDD90183
Air inlet filters (pack/4)	1060318

Please visit www.philips.com/respironics



© 2011 Koninklijke Philips Electronics N.V.

All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

MCI 4103908 1080291 / bdc13037 * JAN 2011

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845



Respironics New Jersey, Inc.
5 Wood Hollow Road
Parsippany, NJ 07054 USA

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics