



Actual size: 140mm (Width) x 175mm (Height)

Ruler is for scaling purposes.

Use this fitting template as a guide to size your patients.


How to use the fitting template:

- Cut out the template along the dotted line
- Follow directions on template to determine the recommended cushion size.
- Please ensure the template is printed to scale.

**Quattro™ Air
AirFit™ F10**


Full face mask
Full face mask **for Her**

Rx Only 
0123




Quattro™ Air
AirFit™ F10


Full face mask
Full face mask **for Her**





SIZE

L
M
S
XS

Rx Only 
0123







Sample size check