

Humidifier/device compatibility list

This supplement specifies the devices compatible with ResMed's humidifiers. For full details regarding the correct use of these devices, please refer to the relevant User Guides.

Note: Not all products are available in all regions.

Device	Humidifier		Heated Humidifiers					External Humidifier
	Passover	HumidAire 2iC™	HumidAire 2i™	HumidAire 3i™ H3i Plus™	H4i™	H4i™ for Stellar	H5i™	HumiCare™ D900
AutoSet CS™ 2		✓	✓					
AutoSet™ CS (S9)							✓	
AutoSet CS™ (S9 PaceWave), AutoSet CS-A™ (S9 PaceWave)							✓	✓
AutoSet Spirit™	✓	✓	✓					
Auto 25 (S9)							✓	
S7™ Series	✓	✓	✓					
S8 AutoScore™, S8 AutoSet Vantage™, S8 AutoSet Spirit™				✓	✓			
S8 Escape™, S8 Lightweight™	✓			✓	✓			
S8 Escape™ II				✓	✓			
S8 Elite™ II, S8 AutoSet Spirit™ II				✓	✓			
S9 AutoSet™, S9 Elite™, S9 Escape™, S9 Escape™ Auto							✓	
VPAP™ III Series	✓	✓	✓					
VPAP™ III ST-A with QuickNav		✓	✓					
VPAP™ Adapt SV		✓	✓					
VPAP™ Adapt (S9)							✓	
VPAP™ Adapt (S9 PaceWave)							✓	✓
V-Auto (S9)							✓	
VPAP™ S (S9), VPAP™ ST (S9), VPAP™ ST-A (S9), VPAP™ ST (S9 iVAPS), VPAP™ ST-A (S9 iVAPS)							✓	✓
VPAP™ IV					✓			
VPAP™ Tx			✓					
VPAP Malibu™			✓					
Stellar™						✓		✓
VS III™								✓
Elisée™								✓

Product	Manufacturer	Distributed by
Stellar	ResMed Germany Inc. Fraunhoferstr. 16 82152 Martinsried Germany	ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA, ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia
VS III, Elisée	ResMed Paris 240 rue de la Motte 77550 Moissy-Cramayel, France	ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia
All others	ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia	ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA, [EC] ResMed (UK) (for Ltd) 96 Jubilee Ave Milton Park Abingdon Oxfordshire OX14 4RW UK

See www.resmed.com for other ResMed locations worldwide. S8, S9, Stellar, PaceWave and VPAP are trademarks of ResMed Ltd and VPAP is registered in U.S. Patent and Trademark Office. Elisée and VS III are trademarks of ResMed Paris SAS and Elisée is registered in U.S. Patent and Trademark Office.
© 2013 ResMed Ltd. 1012241-Eng/9 2013-06